





















## MINI MAINTENANCE KIT

No.	Product Name	Unit	Qty'	Description / Specifications	Pictures
A	<b>Oil Replacement Kit</b>	Box	1	Regular Maintenance Kit (Required a replacement in 12 months or 1250 cycles which comes first)	
A-1	Oil for Replacement	EA	1	Oil for Replacement LB100 (280ml)	
A-2	Waste Oil Bottle	EA	1	Waste Oil Bottle 300ml	
A-3	Oil Injection Syringe	EA	1	Oil Injection Syringe 1EA	
A-4	Oil Injection Tube (Transparent)	EA	1	Oil Injection Tube (Transparent) Translucent, Φ8, 220mm	
A-5	Oil Injection Adaptor (White)	EA	1	Oil Injection Adaptor (Transparent) Female Luer lock connector, PP, 4.80mm	
A-6	Oil Extraction Syringe	EA	1	Oil Extraction Syringe 1EA	
A-7	Oil Extraction Tube (Black)	EA	1	Oil Extraction Tube (Black) Black, φ6, 220mm	
A-8	Oil Extraction Adaptor (Black)	EA	1	Oil Extraction Adaptor (Transparent) Female Luer lock connector, PP, 4.00mm	

No.	Product Name	Unit	Qty'	Description / Specifications	Pictures
B	<b>Filter Replacement Kit (6M)</b>	Box	1	Regular Maintenance Kit (Required a replacement in 6 months or 625 cycles which comes first)	
B-1	HEPA Filter	EA	1	HEPA Filter Storage temperature: + 5~ 40°C	
B-2	Ozone Filter	EA	1	Ozon Filter Shelf life: None/Storage temperature: + 5~ 80°C(Included Ozone Filter Hand Knob 2ea)	
B-3	Chamber Filter	EA	1	Chamber Hole Filter Shelf life: None Storage temperature: + 5~ 80°C	
B-4	VC Cartridge	EA	2	VC Cartridge Shelf life: None / Not included hydrogen peroxide Required a calibration every 6 months or when the installation place is changed	

No.	Product Name	Unit	Qty'	Description / Specifications	Pictures
C	<b>Filter Replacement Kit (12M)</b>	Box	1	Regular Maintenance Kit (Required a replacement in 12 months or 1250 cycles which comes first)	
C-1	HEPA Filter	EA	1	HEPA Filter Storage temperature: + 5~ 40°C	
C-2	Ozone Filter	EA	1	Ozon Filter Shelf life: None/Storage temperature: + 5~ 80°C (Included Ozone Filter Hand Knob 2ea)	
C-3	Chamber Filter	EA	1	Chamber Filter Shelf life: None Storage temperature: + 5~ 80°C	
C-4	VC Cartridge	EA	2	VC Cartridge Shelf life: None / Not included hydrogen peroxide Required a calibration every 6 months or when the installation place is changed	
C-5	Deodorization Filter	EA	1	Deodorization Filter Shelf life: None Storage temperature: + 5~ 80°C (Included O-ring 2ea)	

\* It requires to be replaced every 6 months or 625 cycles which comes first, and 12 months or 1250 cycles which comes first.

\*\* Maintenance Kit is not cover under warranty.